



Recent Advances of Laser based Trace Gas Monitoring Technology: Challenges and Opportunities

Frank K. Tittel

Rice Quantum Institute, Rice University, Houston, TX

<http://ece.rice.edu/lasersci/>

OUTLINE

Physics
Seminar

Lund,
Sweden
Nov. 23,
2004

- Motivation and Technology Issues
- Fundamentals of Laser Absorption Spectroscopy
- Selected Applications of Trace Gas Detection
 - Off Axis-ICOS based detection of NO
 - Direct Absorption Spectroscopy of OCS
 - Quartz Enhanced Laser-PAS of N₂O and H₂CO
 - CO₂ Flux and Isotopic Ratio Measurements
- Conclusions and Outlook

Work supported by NASA, NSF, NIH and Welch Foundation

Acknowledgements

M. Allen (PSI)

F. Capasso (Harvard)

J. Faist (Neuchatel)

A. Fried (NCAR)

C. Gmachl (Princeton)

J.C. Graf (NASA-JSC)

R. Yang (NASA-JPL)

P. Frantz (Rice)

M. Fraser (Rice)

D. Richter (NCAR)

- National Aeronautics and Space Administration
- National Institute of Health
- National Science Foundation
- Texas Advanced Technology Program
- Welch Foundation



1996 Nobel Laureates in Chemistry, R.F. Curl, H. Kroto & R.E. Smalley



Motivation: Wide Range of Gas Sensing Applications


- **Urban and Industrial Emission Measurements**
 - Industrial Plants
 - Combustion Sources and Processes (eg. early fire sensing)
 - Automobile and Aircraft Emissions
- **Rural Emission Measurements**
 - Agriculture and Animal Facilities
- **Environmental Monitoring**
 - Atmospheric Chemistry
 - Volcanic Emissions
- **Chemical Analysis and Industrial Process Control**
 - Chemical, Pharmaceutical, Food & Semiconductor Industry
- **Spacecraft and Planetary Surface Monitoring**
 - Crew Health Maintenance & Human Life Support Program
- **Medical Diagnostics** (eg. breath analysis)
- **Biohazard and Toxic Chemical Detection**
- **Fundamental Science, Photochemistry and Life Sciences**




International Space Station



Existing Methods for Trace Gas Detection



Direct Laser Absorption Spectroscopy




Beer-Lambert's Law of Linear Absorption

$$I(\nu) = I_0 \cdot e^{-\alpha(\nu) \cdot P_a \cdot L}$$

$\alpha(\nu)$ - absorption coefficient [$\text{cm}^{-1} \text{atm}^{-1}$]; L - path length [cm]

ν - frequency [cm^{-1}]; P_a - partial pressure [atm]



$$\alpha(\nu) = C \cdot S(T) \cdot g(\nu - \nu_0)$$

C - total number of molecules of absorbing gas/atm/ cm^3 [$\text{molecule} \cdot \text{cm}^{-3} \cdot \text{atm}^{-1}$]

S - molecular line intensity [$\text{cm} \cdot \text{molecule}^{-1}$]

$g(\nu - \nu_0)$ - normalized lineshape function [cm], (Gaussian, Lorentzian, Voigt)



Sensitivity Enhancement Techniques

- **Optimum Molecular Absorbing Transition**
 - Overtone or Combination Bands (NIR)
 - Fundamental Absorption Bands (MID-IR)
- **Long Optical Pathlength**
 - Multipass Absorption Cell (White, Herriot)
 - Cavity Enhanced and Cavity Ringdown Spectroscopy
 - Open Path Monitoring (with retro-reflector)
 - Fiberoptic Evanescent Wave Spectroscopy
- **Spectroscopic Detection Schemes**
 - Frequency or Wavelength Modulation
 - Balanced Detection
 - Zero-air Subtraction
 - Photoacoustic Spectroscopy
 - Noise Immune Cavity Enhanced-Optical Heterodyne Molecular Spectroscopy (NICE-OHMS)



CW IR Source Requirements for Laser Spectroscopy

REQUIREMENTS


- Sensitivity (% to ppt)
- Selectivity
- Multi-gas Components
- Directionality
- Rapid Data Acquisition
- Room Temperature
- Field deployable

IR SOURCE



- Power
- Narrow Linewidth
- Tunable Wavelengths
- Beam Quality
- Fast Time Response
- No Consumables
- Compact & Robust

9


IR Laser Sources and Wavelength Coverage



HITRAN Simulation of Absorption Spectra (3.1-5.5 & 7.6-12.5 μm)



State of the art: QCL performances



Kindly provided by Prof. J. Faist, University of Neuchatel, Nov. 2004

Key Characteristics of Mid-IR Quantum Cascade Lasers for Spectroscopy

- QC laser wavelengths cover entire mid-IR range from 3.3 to 24 μm determined by thickness of the quantum well and barrier layers of the active region
- Intrinsically high power lasers (determined by number of stages of injector-active quantum well gain regions)
 - CW: ~ 100 mW @ 80°K and 1- 640 mWs @ 295 °K
 - Pulsed: >1 W peak at room temperature, ~ 50 mW avg. @ 0 °C (up to 80 % duty cycle) to 100 mWs (56% d.c)
- High spectral purity: $< \text{kHz} - 330\text{MHz}$ ($.01$ to $.001$ cm^{-1})
- Wavelength tunable by current (~ 1 cm^{-1}) or temperature scanning (~ 10 cm^{-1}); ~ 160 cm^{-1} with external cavity grating
- High reliability: long lifetime, robust operation and reproducible emission wavelengths




11


Important Biomedical Target Gases

Molecule	Formula	Biological/Pathology Indication
Pentane	$\text{CH}_3(\text{CH}_2)_3\text{CH}_3$	Lipid peroxidation, oxidative stress associated with inflammatory diseases, transplant rejection, breast and lung cancer
Ethane	C_2H_6	Lipid peroxidation and oxidative stress
CO_2 isotope ratio	$^{13}\text{CO}_2 / ^{12}\text{CO}_2$	Marker for Heliobacter pylori infection, Gastrointestinal and hepatic function
Carbonyl Sulfide	COS	Liver disease and acute rejection in lung transplant recipients (10-500 ppb?)
Carbon disulfide	CS_2	Schizophrenia
Ammonia	NH_3	Hepatic encephalopathy, liver and renal diseases, fasting response
Formaldehyde	HCHO	Cancerous tumors, breast cancer (400-1500 ppb)
Nitric Oxide	NO	Inflammatory and immune responses (e.g., asthma) and vascular smooth muscle response (6-100 ppb)
Hydrogen Peroxide	H_2O_2	Airway Inflammation, Oxidative stress (1-5 ppb)
Carbon Monoxide	CO	Smoking response, CO poisoning, vascular smooth muscle response, platelet aggregation (400-3000 ppb)
Ethylene	$\text{H}_2\text{C}=\text{CH}_2$	Oxidative stress, cancer
Acetone	CH_3COCH_3	Fasting response, diabetes mellitus response, ketosis

Nitric Oxide Concentrations in Exhaled Breath by Healthy Humans





Off-Axis Integrated Cavity Output Spectroscopy Based Gas Sensor



- * Novel compact gas cell design of length: 3.8 – 5.3 cm and cell volumes < 80 cm³;
- * Low loss mirrors (ROC 1m): ~ 60 -250 ppm, R=99.975, L_{eff} =170-800 m
- * Rapid eNO concentration measurements during a single breath cycle are feasible

Off-axis ICOS Detection of NO



- 95 and 490 ppb NO/N₂ calibration mixture at 100 Torr total pressure
- Effective optical path ~ 70 m (1,350 passes)

Voigt fit of measured NO absorption line at 1920.7 cm⁻¹ for a concentration of 95 ppb

Noise-equivalent sensitivity is 10 ppb for 1σ deviation of the best fit coefficient.
 Detection sensitivity: $1.0 \times 10^{-7} \text{ cm}^{-1} \text{ Hz}^{-1/2}$




17

NO from Nasal Exhaled Air (OA-ICOS and wavelength modulation spectroscopy)



Medically approved collection bag for human breath samples



Averaged 2f signal of the OA-ICOS cavity output


- (1) Nasal NO concentration: 53 ppb
- (2) 95 ppb NO/N₂ calibration mixture



Exhaled Carbonyl Sulfide



- A 2001 study by the T. H. Risby group at John Hopkins University demonstrated that elevated levels of COS could have a diagnostic role in the detection of acute allograft rejection in lung transplant recipients
 - S. M. Studer, et. al., J. of Heart and Lung Transplantation, 20(11), 1158-66 (2001).

Measured with gas chromatography and flame ionization detection




17

OCS Sensor Architecture






G. Wysocki et al. Applied Optics 43, 8040-8046, 2004

OCS ro-vibrational ν_3 spectrum at $\sim 4.85\mu\text{m}$




OCS Concentration Calibration of QCL Sensor








Schematic of a Portable LAS System




QC laser based measurements of CO trace gas above cell cultures



- Measured CO production rates of viable cultures of vascular smooth muscle cells
- Achieved a detection limit of for CO of ~ 20 ppb


Comparative Size of Absorbance Detection Modules (ADM)



Merits of QE Laser-PAS based Trace Gas Detection

- Immune to ambient and flow acoustic noise, laser noise and etalon effects
- Dramatic reduction of sample volume (< 1 mm³)
- High sensitivity (ppm to ppb gas concentration levels) and excellent dynamic range
- Applicable over a wide range of pressures
- Temperature, pressure and humidity insensitive
- Rugged and low cost compared to LAS that requires a multipass absorption cell and infrared detector(s)
- Potential for optically multiplexed concentration measurements


Resonant photoacoustic spectroscopy



27



Quartz-Enhanced Photoacoustic Spectroscopy (QEPAS)




28

A. Kosterev et al. Optics Letters 27, 1902, 2002

QEPAS versus Traditional PAS

Parameter	Traditional PAS	QEPAS
f , Hz	100 to 4000	Presently ~32 760
Q	20 to 200	10 000 to 30 000
Q vs. pressure	INCREASES (high spectral resolution is problematic)	DECREASES (high spectral resolution is achievable)
Sample volume	$>10 \text{ cm}^3$	$<1 \text{ mm}^3$
Sensitivity to ambient acoustic and flow noise	Usually high	None observed
Pathlength involved	~10 cm	(a) 0.3mm, (b) 5mm


Diode Laser based QEPAS Gas Sensor Architecture



32




Infrared NH₃ Absorption Spectra




M. Webber et al. 2003, Pranalytica

Ammonia Detection using a 1.53 μm Telecom Diode Laser




Spectral simulations based on data from Webber et al., APPLIED OPTICS **40**, 2031-2042 (2001)]



QEPAS spectra at different pressures of NH₃-N₂ gas mixture; t=0.3s, 38 mW diode laser excitation power at 6529 cm⁻¹

To appear in Applied Optics (Nov 2004)

Calibration and Linearity of QEPAS based NH₃ Sensor




Noise-equivalent (1s) concentration (NEC). for $\tau=1s$ time constant is 0.65 ppm for 38 mW excitation power

90 last points of each step averaged

Noise-equivalent absorption (NEA) coefficient $k=7.2 \times 10^{-9} \text{ cm}^{-1} \text{ W/Hz}^{1/2}$

Equivalent circuit of a quartz TF



$$\omega_0 = \sqrt{\frac{1}{LC}}$$


$$Q = \frac{1}{R} \sqrt{\frac{L}{C}}$$

$$\sqrt{\langle I_N^2 \rangle} = \sqrt{\frac{4k_B T}{R}}$$

36 "QUARTZ CRYSTAL RESONATORS AND OSCILLATORS For Frequency Control and Timing Applications", tutorial by John R. Vig, U.S. Army Communications-Electronics Command (July 2001)



Noise analysis



Amplifier:
 •Fast
 •Low noise
 •High impedance
 •Low 1/f noise

$$S_1 = \sqrt{4k_B T R_g}; \quad R_g = 10 \text{ M}\Omega \Rightarrow S_1 = 4.1 \cdot 10^{-7} \frac{\text{V}}{\sqrt{\text{Hz}}}$$


$$S_2 = \sqrt{\frac{4k_B T}{R}} R_g; \quad R = 100 \text{ k}\Omega \Rightarrow S_2 = 4.1 \cdot 10^{-6} \frac{\text{V}}{\sqrt{\text{Hz}}} \text{ (at 760 Torr)}$$

$$S = \sqrt{S_1^2 + S_2^2} \approx S_2 \text{ (at resonance)} \quad \text{Noise goes up as } \sqrt{Q}.$$

37




Resonant Response, Q=116 600



38





Pressure Dependence of Q Factor of a Typical TF




39



QCL based Quartz-Enhanced Photoacoustic Spectrometer



N₂O Detection in Ambient Air at 4.6μm (2195.6 cm⁻¹)



Noise-equivalent absorption coefficient $k=1.5 \times 10^{-8} \text{ cm}^{-1} \text{ W/Hz}^{1/2}$

To appear in Applied Physics B 2004




Motivation for Precision Monitoring of H₂CO


- Precursor to atmospheric O₃ production
- Pollutant due to incomplete fuel combustion processes
- Potential trace contaminant in industrial manufactured products
- Medically important gas



Mid-IR DFG Based H₂CO Sensor




H₂CO Detection in Ambient Air at 3.53 μm



D. Rehle et al. Applied Physics B72, 947-952 (2001)






H₂CO and O₃ Concentrations at Deer Park, TX for July 20-31, 2003



43



Advanced DFG System for H₂CO Detection



Laser₁ = 1560 nm
 Laser₂ = 1083 nm

$\lambda_{\text{idler}} = 2831.6417 \text{ cm}^{-1}$

A. Fried et al, Development of DFG source in progress

Motivation for a NCAR Biocomplexity Project (2003-2005)

- Study of the Carbon Cycle: Land, Ocean, Atmosphere
- Combustion of fossil fuels increase atmospheric CO₂ levels and in turn global temperatures
- Currently, oceanic and terrestrial biospheres absorb 50% of anthropogenic atmospheric carbon. Will the sinks cease to take up CO₂ ?
- Major sources and sinks have been identified, but their specific mechanisms are still unclear
- Flexible and real time sensor is needed.
- Bring the instrument to the sample: Land, sea, troposphere, stratosphere, planets !



NCAR C-130

Accurate measurements of ^{12/13}CO₂ ratios and other trace gases are key to attain detailed knowledge of the carbon cycle processes

Intercontinental Chemical Transport Experiment 2004


Each CH₂O Point
Replicate Precision
Over 1 min, N = 60

Avg. Channel Std. Dev

$$\sigma_{1\text{min}} = \frac{\sigma_{1\text{sec}}}{\sqrt{n}} \cong \frac{200 \text{ ppt}}{\sqrt{60}} \cong 26 \text{ ppt}$$

Each Point
1-second data points


Concentration determined
by least squares fit of
Calibration spectrum
to Sample spectra



INTEX July 2004, data provided by Alan Fried et. al., NCAR
Boulder, CO

29


HITRAN Based Simulation of a H₂CO-H₂O-CH₄ Spectrum in Tuning Range of a 3.53μm IC Laser



- H₂CO : 10 ppb
- H₂O : 3%
- CH₄ : 2 ppm
- Optical path: 100 m
- Total pressure: 30 Torr



2f - QEPAS based H₂CO signal at 3.53 μm (2832.48 cm⁻¹)




- [H₂CO]: 13.27 ppm
- Sensitivity: $2.2 \times 10^{-8} \text{ cm}^{-1} \text{ W} / \sqrt{\text{Hz}}$
- QEPAS NE sensitivity for NH₃: $7.2 \times 10^{-9} \text{ cm}^{-1} \text{ W} / \sqrt{\text{Hz}}$
(For comparison)

Applied Physics B 79, 799-803, 2004



IC Laser based Formaldehyde Calibration Measurements with a Gas Standard Generator



- H₂CO absorption frequency: 2832.5 cm⁻¹
- Lock-In time constant: 10 s
- QEPAS parameters
 - Resonance frequency: 32.760 KHz
 - Q-factor: 17336
 - Pressure: 200 Torr
 - Gas Flow: 75 sccm
 - IC laser power: 6 mW



QEPAS Performance for 7 Trace Gas Species (2004)

Molecule (Host)	Frequency, cm ⁻¹	Pressure, Torr	NNEA, cm ⁻¹ W/Hz ^{1/2}	Power, mW	NEC (τ=1s), ppmv
NH ₃ (N ₂)	6528.76	60	7.2×10 ⁻⁹	38	0.65
H ₂ O (exhaled air)	6541.29	90	8×10 ⁻⁹	5.2	580
CO ₂ (exhaled air)	6514.25	90	1.0×10 ⁻⁸	5.2	890
N ₂ O (air+5%SF ₆)	2195.63	50	1.5×10 ⁻⁸	19	0.007
CO (N ₂)	2196.66	50	5.3×10 ⁻⁷	13	0.5
CO (propylene)	2196.66	50	7.4×10 ⁻⁸	6.5	0.14
CH ₂ O (air)	2832.48	200	2.2×10 ⁻⁸	3.4	0.55

NNEA – normalized noise equivalent absorption coefficient.

NEC – noise equivalent concentration for available laser power and τ=1s time constant.



Motivation for Measuring $^{13}\text{CO}_2/^{12}\text{CO}_2$ Isotopic Ratios


- Atmospheric Chemistry: Environmental monitoring of C_y gases ($\text{CO}_2, \text{H}_2\text{O}, \text{CO}, \text{N}_2\text{O}, \text{CH}_4$)
 - Global warming studies
 - Temporal and spatial variations of isotopic ratios
 - Identification of carbon sources and sinks
 - Global carbon budget studies
- Study of planetary gases (e.g. for Mars: $\text{CO}, \text{CO}_2, \text{H}_2\text{O}, \text{CH}_4, \text{O}_3, \text{OCS}$)
- Volcano eruption forecasting and gas emission studies ($\text{CO}_2, \text{HCl}, \text{SO}_2, \text{HF}, \text{H}_2\text{S}, \text{CO}, \text{H}_2\text{O}$)
- Geochemistry
- Medical applications (non-invasive human health monitoring)



51

Volcanological applications

- CO_2 the most abundant component of volcanic gases after H_2O
- $\delta^{13}\text{C}$ is a sensitive tracer of magmatic vs. hydrothermal or groundwater contributions to volcanic gases
- Monitoring $\delta^{13}\text{C}$ can be used in eruption forecasting and volcanic hazard assessment



Isotopic Ratio Measurement Techniques

- Current standard technique: Isotope Ratio Mass Spectrometry (IRMS) $\Delta\delta\sim 0.01-0.05\text{ ‰}$
 - Small mass differences are difficult to measure
 - Not real time
 - Not field deployable
 - Complex sample preparation and sample destruction
- Fourier Transform Spectroscopy $\Delta\delta\sim 0.1-0.2\text{ ‰}$
 - Not selective for compact and intermediate sized platform
- Tunable Laser Absorption Spectroscopy $\Delta\delta\sim 0.2\text{ ‰}$
 - Lead salt lasers
 - Difference Frequency Generation
 - Near-infrared diode
 - Mid-infrared quantum cascade lasers

CO₂ Absorption Line Selection Criteria

- Three strategies:
 - Similar strong absorption of ¹²CO₂ and ¹³CO₂ lines
 - Very sensitive to temperature variations
 - Similar transition lower energies
 - Requires a dual path length approach to compensate for the large difference in concentration between major and minor isotopic species-or-
 - Can be realized if different vibrational transitions are selected for the two isotopes (4.35 μm for ¹³CO₂ and 2.76 μm for ¹²CO₂)*
- For the first 2 strategies both absorption lines must lie in a laser frequency scan window
- Avoid presence of other interfering atmospheric trace gas species



* Proposed scheme by Curl, Uehara, Kosterev and Tittel, Oct. 2002




Comparison of CO₂ line selection and strategy for different current US mid-IR laser-based isotopic ratiometers

Group	Technology	Frequency 12/13 [cm ⁻¹]	δT [K]	Precision
NCAR, UC and Rice U. A. Fried et al; Erdelyi	DFG with NIR TDLs and fiber amplifiers	2299.642 2299.795	0.005	0.8 ‰*
Aerodyne, Harvard U. M. Zahniser et al.	Direct Scan PbSalt TDL, QCL, DFG; Dual optical paths	2314.304 2314.408	0.213	0.2 ‰
Physical Sciences D. Sonnenfroh et al	QCL	2318.1		0.5 to 1‰
Rice University Tittel et al	QCL Dual optical paths	2311.105 2311.399	181 Very large	<1 ‰
U. of Utah Bowling, Picarro	PbSalt TDLs Campbell Scientific Instrum.	2308.225 2308.171	0.006	0.2 ‰
JPL C. Webster	TDLs and QCL, LAS	2303.7 2303.5	0.007	TBD ‰
NASA-Ames Becker et al; Jost, LGR	Direct Scan PbSalt TDLs & QCLs with ICOS	2291.542 2291.680	0.004	4 ‰


HITRAN based Simulation of Ro-vibrational bands suitable for ¹²CO₂/¹³CO₂ ratio measurements




Mid-IR DFG Isotopic Ratio Spectrometer



Current Status of DFG based IRS System at NCAR



QC laser based Isotopic Ratio Sensor at Rice




Bread board: 12x18" (30x45 cm)

The sensor must be operated in a dry nitrogen atmosphere to eliminate atmospheric CO₂ background


59



Dual path length gas cell design for infrared ratio spectrometry




$^{16}\text{O}^{12}\text{C}^{18}\text{O}$ Spectra at 2320.2 cm^{-1}



D. Weidmann et al. Applied Physics B 2004



DFG of $^{13}\text{CH}_4/^{12}\text{CH}_4$ at $3.3\text{ }\mu\text{m}$




ADVANTAGES:

WIDELY TUNABLE
NARROW LINEWIDTH $<0.001\text{ cm}^{-1}$
EXCELLENT BEAM QUALITY

DISADVANTAGE:

COMPLEXITY

AERODYNE, Inc. M. Zahniser et. al.
SANDIA NATIONAL LABORATORY




Towards CW Operation of DFB QCLs

- Main drawbacks of QCL pulsed operation
 - Pulse to pulse intensity variation
 - Linewidth broadening by thermal chirping
 - Requirement of nanosecond electronics
- Efforts towards achieving quasi-RT CW DFB QCLs
 - M. Beck *et al.*, *Science*, **295**, 301-305, 2002
 - T. Aellen *et al.*, *Applied Physics Letters*, **83**, 1929, 2003
 - A. Evans, *et al.*, *Applied Physics Letters*, **84**, 314, 2004 [FP QCL]
 - M. Troccoli *et al.*, 5th Workshop on Quantum Cascade Lasers, Freiburg, Sept. 23-24 2004





CW RT DFB QC laser based NH₃ Sensor @ 9 μm (1113 cm⁻¹)



Wavelength Modulation Spectroscopy of NH_3

- QCL Drive Current :
Quasi CW +
Wavelength modulation





- Calibration with a
1038 ppm $\text{NH}_3:\text{N}_2$ mixture
 1σ extrapolated sensitivity
82 ppb.m/ $\sqrt{\text{Hz}}$

⇒ Improvement by a factor of
3 compared to direct
absorption spectroscopy

67

λ -Tunable Pulsed QC Laser GCEC Architecture



H. Le et al., Univ. of Houston

NASA Atmospheric & Mars Gas Sensor Platforms



Tunable laser sensor for earth's stratosphere

Aircraft laser absorption spectrometer



Tunable laser planetary spectrometer



Conclusions and Future Directions

- **Quantum Cascade Laser based Trace Gas Sensors**
 - Compact and robust sensors based on QC-LAS and QE L-PAS
 - High sensitivity (10^{-4} - 10^{-5}) and selectivity (3 to 500 MHz)
 - Dramatic reduction of sample volume ($\sim 0.2 \text{ mm}^3$)
 - Detected trace gases: NH_3 , CH_4 , N_2O , CO_2 , CO , NO , H_2O , COS , C_2H_4 , $\text{C}_2\text{H}_5\text{OH}$, SO_2 , H_2CO and several isotopic species of C, O, N and H.
- **Applications in Trace Gas Detection**
 - Environmental monitoring (NH_3 , CO , CH_4 , C_2H_4 , N_2O , CO_2)
 - Industrial process control and chemical analysis (NO , NH_3)
 - Medical Diagnostics (NO , CO , COS , CO_2 , C_2H_4)
- **Future Directions and Collaborations**
 - Cavity enhanced (ICOS) and QE L-PAS spectroscopy based applications using novel thermoelectrically cooled cw and broadly wavelength tunable quantum cascade lasers
 - Applications using new near IR interband and far-IR intersub-band quantum cascade lasers